


Zebrafish Prpf18 Antibody / Pre-mRNA-splicing factor 18 (RZ1281)

Catalog No.	Formulation	Size
RZ1281	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q6GMH0
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Prpf18 antibody is available for research use only.



Western blot analysis of Prpf18 protein using Zebrafish Prpf18 antibody and 1) whole female zebrafish tissue lysates and 2) whole male zebrafish tissue lysates. Predicted molecular weight ~39 kDa.