


Zebrafish Prpf6 Antibody / Pre-mRNA processing factor 6 (RZ1283)


Catalog No.	Formulation	Size
RZ1283	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	A0A2R8PWG5
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunofluorescence : 5ug/ml
Limitations	This Zebrafish Prpf6 antibody is available for research use only.



Immunofluorescent staining of PRPF6 protein using Zebrafish PRPF6 antibody (green).
HIER: steam section in pH8 EDTA buffer for 20 min.



Western blot analysis of Prpf6 protein using Zebrafish Prpf6 antibody and zebrafish embryo tissue lysates. Predicted molecular weight ~109 kDa.

