


Zebrafish Psma3 Antibody / Proteasome subunit alpha type 3 (RZ1285)

Catalog No.	Formulation	Size
RZ1285	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q4V918
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Psma3 antibody is available for research use only.



Western blot analysis of Psma3 protein using Zebrafish Psma3 antibody and 1) zebrafish head, 2) whole female zebrafish and 3) whole male zebrafish tissue lysate. Predicted molecular weight ~28 kDa.