


Zebrafish Psmd7 Antibody / 26S proteasome non-ATPase regulatory subunit 7 (RZ1295)

Catalog No.	Formulation	Size
RZ1295	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q7ZYX7
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Psmd7 antibody is available for research use only.



Western blot analysis of Psmd7 protein using Zebrafish Psmd7 antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~37 kDa.