


Zebrafish Rbbp4 Antibody / Retinoblastoma-binding protein 4 (RZ1301)

Catalog No.	Formulation	Size
RZ1301	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q6P3H7
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Rbbp4 antibody is available for research use only.



Western blot analysis of Rbbp4 protein using Zebrafish Rbbp4 antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~48 kDa, commonly observed at 48-55 kDa. (human similarity)