


Zebrafish Sod2 Antibody / Superoxide dismutase 2 (RZ1312)

Catalog No.	Formulation	Size
RZ1312	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Predicted Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q6P980
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Sod2 antibody is available for research use only.



Western blot analysis of Sod2 protein using Zebrafish Sod2 antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~25 kDa.