


Zebrafish Sox3 Antibody / SRY-box 3 (RZ1194)

Catalog No.	Formulation	Size
RZ1194	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q6EJB7
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Sox3 antibody is available for research use only.



Western blot analysis of Sox3 protein using Zebrafish Sox3 antibody and 1) whole female zebrafish tissue lysates and 2) zebrafish embryo tissue lysates. Predicted molecular weight ~33 kDa.