


Zebrafish Spi1b Antibody | Pu.1b | Spi1 (RZ1314)

Catalog No.	Formulation	Size
RZ1314	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	A0A8M2BG46
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Spi1b antibody is available for research use only.



Western blot analysis of Spi1b protein using Zebrafish Spi1b antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~35 kDa.