


Zebrafish Star Antibody / Steroidogenic acute regulatory protein (RZ1317)

Catalog No.	Formulation	Size
RZ1317	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q9DG10
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Star antibody is available for research use only.



Western blot analysis of Star protein using Zebrafish Star antibody and zebrafish head tissue lysates. Predicted molecular weight ~32 kDa. The human preprotein is ~37 kDa with a 30-32 kDa mature form.