


Zebrafish Tbx5a Antibody / Tbx5 / T-box protein 5a (RZ1320)

Catalog No.	Formulation	Size
RZ1320	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q9IAK8
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Tbx5a antibody is available for research use only.



Western blot analysis of Tbx5a protein using Zebrafish Tbx5a antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~55 kDa.