


Zebrafish Tfap2a Antibody / Ap2 alpha (RZ1190)

Catalog No.	Formulation	Size
RZ1190	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	A0A2R8QRQ5
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Tfap2a antibody is available for research use only.



Western blot analysis of Tfap2a protein using Zebrafish Tfap2a antibody and 1) zebrafish head tissue lysates and 2) zebrafish embryo tissue lysates. Predicted molecular weight ~50 kDa.