


Zebrafish Tnfa Antibody / Tnf alpha (RZ1321)

Catalog No.	Formulation	Size
RZ1321	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q08CQ3
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Tnfa antibody is available for research use only.



Western blot analysis of Tnfa protein using Zebrafish Tnfa antibody and 1) zebrafish head, 2) whole female zebrafish and 3) whole male zebrafish tissue lysate. Predicted molecular weight ~26 kDa.