


ZFP64 Antibody (F40052)

Catalog No.	Formulation	Size
F40052-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40052-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9NTW7
Applications	Western Blot : 1:1000
Limitations	This ZFP64 antibody is available for research use only.



Description

ZFP64 may be involved in transcriptional regulation.

Application Notes

Titration of the ZFP64 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 525-551 from the human protein was used as the immunogen for this ZFP64 antibody.

Storage

Aliquot the ZFP64 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.