


ZMPSTE24 Antibody (R33968)

Catalog No.	Formulation	Size
R33968-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10269
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ZMPSTE24 antibody is available for research use only.



Western blot testing of human skin lysate with ZMPSTE24 antibody at 0.5ug/ml.
Predicted molecular weight: ~55kDa.

Description

Additional name(s) for this target protein: Zinc metalloproteinase (STE24 homolog, yeast); FACE1

Application Notes

Optimal dilution of the ZMPSTE24 antibody should be determined by the researcher.

Immunogen

Amino acids ERLQALKTMKQH were used as the immunogen for this ZMPSTE24 antibody.

Storage

Aliquot and store the ZMPSTE24 antibody at -20oC.