


ZNF322B Antibody (F40011)

Catalog No.	Formulation	Size
F40011-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40011-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q6U7Q0
Applications	Western Blot : 1:1000
Limitations	This ZNF322B antibody is available for research use only.



Description

May be involved in transcriptional regulation (By similarity).

Application Notes

Titration of the ZNF322B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 326-354 from the human protein was used as the immunogen for this ZNF322B antibody.

Storage

Aliquot the ZNF322B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.